



# Super Edge SNAP<sup>®</sup> Laboratory Glass Columns

*“Next Generation” technology for high-performance preparative chromatography*



Colum ID mm	Length mm	Pressure MPa	Version	Aqueous buffer (AB)	Solvent Resistant (SR)
10	125	4.0	Body	Borosilicate glass	
15	250	3.5	Piston	Acetal	PEEK
25	350	2.4	Frit	Polyethylene	Stainless Steel
35	500	1.8	O-Ring	EPDM	Viton or Kalrez <sup>®</sup>
50	750	1.3	Temp. range	4 - 40 deg. C	4 - 80 deg. C
	1000				

**JNC JNC CORPORATION**

JNC Corporation Life Chemicals Launch Office  
2-1, Otemachi 2-Chome, Chiyoda-ku,  
Tokyo 100-8105 Japan

SNAP<sup>®</sup> is a registered US Trademark of ELS

**PARTNERS**



**MR SANGHAVI & CO.**

Trivedi Bldg., 2 Ansari Rd., Vile Parle(W), Mumbai 400056  
Tel +91-250-2390409 ; Fax +91-250-2390112  
Sanghavi.rajiv@gmail.com